


Product Datasheet

GARP antibody (orb222768)

Description	Mouse monoclonal antibody to GARP
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Purified human sGARP protein
Target	LRRC32
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Western blotting: Recommended dilution: 1-2 µg/ml; positive control: human thrombocytes, reducing conditions
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	Q14392
Entrez	2615
Dilution Range	Western blotting: Recommended dilution: 1-2 µg/ml; positive control: human thrombocytes, reducing conditions
Expiration Date	12 months from date of receipt.



Western blotting analysis (reducing cond...