



Product Datasheet

CD42a antibody (orb222771)

Description	Mouse monoclonal antibody to CD42a
Reactivity	Canine, Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Human acute lymphoblastic leukemia cells
Target	CD42a
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P14770
Entrez	2815
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of human thrombocytes (red-fi...



Flow cytometry surface staining pattern ...