


Product Datasheet


CD266 antibody (orb222780)

Description	Mouse monoclonal antibody to CD266
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IHC-Fr, WB
Immunogen	human CD266-transfected P815 cells
Target	CD266
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml
Isotype	Mouse IgG2b
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	Q9NP84
Entrez	51330
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml
Expiration Date	12 months from date of receipt.



Separation of HUVEC cells stained using

...



Flow cytometry surface staining pattern ...