


Product Datasheet


CD266 antibody (orb222781)

Description	Mouse monoclonal antibody to CD266
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FA, FC, IHC-Fr, WB
Immunogen	human CD266-transfected P815 cells
Target	CD266
Preservatives	Phosphate buffered saline (PBS), pH 7.4
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Low Endotoxin, Azide free format. Product does not contain preservatives and is intended primarily for functional applications. Preservative containing buffers are available for IHC applications
Isotype	Mouse IgG2b
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	Q9NP84
Entrez	51330
Dilution Range	Functional application: Blocking. Flow cytometry: Recommended dilution: 1-4 µg/ml
Expiration Date	12 months from date of receipt.



Separation of HUVEC cells stained using

...



Flow cytometry surface staining pattern ...