



Product Datasheet

CD235a antibody (orb222783)

Description	Mouse monoclonal antibody to CD235a
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC-Fr, IHC-P
Immunogen	Membrane preparation from splenic hairy cell leukemia
Target	CD235a
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Immunohistochemistry (paraffin sections): Recommended dilution: 10 µg/ml. Flow cytometry: Recommended dilution: 1-4 µg/ml
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P02724
Entrez	2993
Dilution Range	Immunohistochemistry (paraffin sections): Recommended dilution: 10 µg/ml. Flow cytometry: Recommended dilution: 1-4 µg/ml
Expiration Date	12 months from date of receipt.



Flow cytometry analysis (surface staining)



Immunohistochemistry staining of human aorta