



Product Datasheet

CD11a antibody (orb222795)

Description	Mouse monoclonal antibody to CD11a
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FA, FC, IP
Immunogen	Leukocytes from a patient suffering from a LGL-type leukaemia.
Target	CD11a
Preservatives	Phosphate buffered saline (PBS), pH 7.4
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Low Endotoxin, Azide free format. Product does not contain preservatives and is intended primarily for functional applications. Preservative containing buffers are available for IHC applications
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P20701
Entrez	3683
Dilution Range	Functional application: The antibody MEM-25 partially blocks binding of LFA-1 complex to ICAM-1. Immunoprecipitation: Excellent antibody for immunoaffinity purification of LFA-1 complex. Flow cytometry: Recommended dilution: 2 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of human monocytes (red-filled)



Flow cytometry surface staining pattern ...