



Product Datasheet

CD41 antibody (orb222812)

Description	Rat monoclonal antibody to CD41
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC-Fr, IHC-P, IP
Immunogen	Murine platelets
Target	CD41
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Isotype	Rat IgG1 kappa
Clonality	Monoclonal
Purity	Purified by protein-G affinity chromatography.
Uniprot ID	Q9QUM0
Entrez	16399
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of murine CD41 positive throm...



Flow cytometry surface staining pattern ...