



Product Datasheet

CD178 antibody (orb222814)

Description	Mouse monoclonal antibody to CD178
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IP, WB
Immunogen	L5178Y mouse T lymphoma cells expressing recombinant human CD178
Target	CD178
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P48023
Entrez	356
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml
Expiration Date	12 months from date of receipt.



Separation of FasL transfected L5178Y ce...



Flow cytometry surface staining pattern ...