


Product Datasheet

CD157 antibody (orb222821)

Description	Mouse monoclonal antibody to CD157
Reactivity	Human, Primate
Conjugation	Unconjugated
Tested Applications	FC, IHC-Fr, IP, WB
Immunogen	Human CD157
Target	CD157
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Western blotting: Non-reducing conditions. Immunohistochemistry: Recommended dilution: 5-10 µg/ml. Flow cytometry: Recommended dilution: 1-12 µg/ml
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	Q10588
Entrez	683
Dilution Range	Western blotting: Non-reducing conditions. Immunohistochemistry: Recommended dilution: 5-10 µg/ml. Flow cytometry: Recommended dilution: 1-12 µg/ml
Expiration Date	12 months from date of receipt.



Flow
cytometry
analysis
(surface
stainin...