



Product Datasheet

CD3 antibody (Biotin) (orb222827)

Description	Mouse monoclonal antibody against CD3 conjugated to Biotin
Reactivity	Human, Primate
Conjugation	Biotin
Tested Applications	FC, IHC-Fr, IP
Immunogen	human thymocytes followed by Sezary T cells
Target	CD3
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Extracellular and intracellular staining; recommended dilution: 1-12 µg/ml.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified antibody is conjugated with biotin LC-NHS ester under optimum conditions and unconjugated antibody and free biotin are removed by size-exclusion chromatography.
Uniprot ID	P07766
Entrez	916
Dilution Range	Flow cytometry: Extracellular and intracellular staining; recommended dilution: 1-12 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of human CD3 positive lymphoc...



Flow cytometry surface staining pattern ...