


Product Datasheet

Human Transthyretin protein (orb222833)

Description	Human Transthyretin protein
Conjugation	Unconjugated
Tested Applications	SDS-PAGE, WB
Target	Transthyretin
Form/Appearance	0.14mol/L NaCl, 0.027mol/L KCl, 0.01mol/L Na ₂ PO ₄ , 0.018mol/L KH ₂ PO ₄ , pH 7.4
Concentration	1.0 mg/mL
Storage	Store at 4°C for short terms and -20°C for long terms storage
Note	For research use only
Protein Sequence	MASHRLLLLC LAGLVFVSEA GPTGTGESKC PLMVKVLDVAV RGSPAINVAV HVFRKAADD T WEPFASGKTS ESGELHGLTT EEEFVEGIYK VEIDTKSYWK ALGISPFHEH AEVVFTANDS GPRRYTIAAL LSPYSYSTTA VVTNPKE
Purity	> 85% as determined by SDS-PAGE
MW	21 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human TTR
protein