



Product Datasheet

EGFR (Phospho-Thr693) antibody (orb222836)


Description	Rabbit polyclonal antibody to EGFR
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P00533
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, IF/ICC: 1:100-1:500, ELISA: 1:40000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HuvEc cells...



Western blot analysis of A431 cells using...



Immunohistochemistry analysis of human breast tissue...