



Product Datasheet

I kappa B-beta (Phospho-Ser23) antibody (orb222838)


Description	Rabbit polyclonal antibody to I kappa B-beta
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	N-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q15653
Dilution Range	WB: 1:500-3000, IHC-P: 1:50-100, IF/ICC: 1:100-1:500, ELISA: 1:40000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...



Western blot analysis of HeLa cells usin...



Immunohistochemical analysis of human lu...