



Product Datasheet

Keratin 8 (Phospho-Ser74) antibody (orb222840)

Description	Rabbit polyclonal antibody to Keratin 8
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P05787
Dilution Range	WB: 1:500-3000, IHC-P: 1:50-100, ELISA: 1:40000
Expiration Date	12 months from date of receipt.



Western blot analysis of 293 cells using...



Immunohistochemical analysis of human br...