



Product Datasheet

Tyrosine Hydroxylase (Phospho-Ser31) antibody (orb222848)

Description	Rabbit polyclonal antibody to Tyrosine Hydroxylase
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	N-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P07101
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, ELISA: 1:1000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of human brain tissue showing Tyrosine Hydroxylase (Phospho-Ser31) staining.



Western blot analysis of 293 cells using Tyrosine Hydroxylase (Phospho-Ser31) antibody.