



Product Datasheet

Tyrosine Hydroxylase (Phospho-Ser40) antibody (orb222849)

Description	Rabbit polyclonal antibody to Tyrosine Hydroxylase
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	N-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P07101
Dilution Range	WB: 1:500-3000, IHC-P: 1:50-100, ELISA: 1:40000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of human brain tissue.



Western blot analysis of RAW264.7 cells.