



Product Datasheet

AKT (Phospho-Thr308) antibody (orb222851)

Description	Rabbit polyclonal antibody to AKT (Phospho-Thr308)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9Y243 , P31751 , P31749
Dilution Range	WB: 1:500-1:3000, ELISA: 1:1000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemistry analysis of human b...



Western blot analysis of COLO using Akt ...