



Product Datasheet

PLCG2 (Phospho-Tyr753) antibody (orb222856)

Description	Rabbit polyclonal antibody to PLCG2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P16885
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, ELISA: 1:5000
Expiration Date	12 months from date of receipt.



Immunohistochemistry analysis of human c...



Western blot analysis of HepG2 cells usi...