



Product Datasheet

CDC25A (Phospho-Ser76) antibody (orb222857)

Description	Rabbit polyclonal antibody to CDC25A
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P30304
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, ELISA: 1:10000
Expiration Date	12 months from date of receipt.



Western blot analysis of HuvEc using CDC...



Immunohistochemistry analysis of human l...