

Product Datasheet

PAK1 (Phospho-Thr212) antibody (orb222874)

Description	Rabbit polyclonal antibody to PAK1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q13153
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, IF/ICC: 1:100-1:500, ELISA: 1:5000
Expiration Date	12 months from date of receipt.