

Product Datasheet

Histone H3 (Phospho-Thr11) antibody (orb222879)

Description	Rabbit polyclonal antibody to Histone H3
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF
Immunogen	N-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q71DI3, P84243, P68431
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, IF/ICC: 1:100-1:500, ELISA: 1:40000
Expiration Date	12 months from date of receipt.