

Product Datasheet

PI3-kinase p85-alpha/gamma (Phospho-Tyr467/199) antibody (orb222945)

Description	Rabbit polyclonal antibody to PI3-kinase p85-alpha/gamma
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P27986, Q92569
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, ELISA: 1:10000
Expiration Date	12 months from date of receipt.