



## Product Datasheet

### Akt2 (Phospho-Ser474) antibody (orb222953)

|                            |                                                     |
|----------------------------|-----------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to Akt2                  |
| <b>Species/Host</b>        | Rabbit                                              |
| <b>Reactivity</b>          | Human, Mouse, Rat                                   |
| <b>Conjugation</b>         | Unconjugated                                        |
| <b>Tested Applications</b> | IHC, WB                                             |
| <b>Immunogen</b>           | C-terminal                                          |
| <b>Concentration</b>       | 1mg/ml                                              |
| <b>Note</b>                | For research use only                               |
| <b>Clonality</b>           | Polyclonal                                          |
| <b>Uniprot ID</b>          | <b>P31751</b>                                       |
| <b>Dilution Range</b>      | WB: 1:500-1:1000, IHC-P: 1:50-1:100, ELISA: 1:20000 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                     |



Western blot analysis of 293 cells using...



Immunohistochemical analysis of human lu...