

Product Datasheet

ATF2 (Phospho-Ser112 or 94) antibody (orb222956)

Description	Rabbit polyclonal antibody to ATF2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P15336
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:100, ELISA: 1:20000
Expiration Date	12 months from date of receipt.