



Product Datasheet

CDC2 (Phospho-Thr161) antibody (orb222968)

Description	Rabbit polyclonal antibody to CDC2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P06493
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, ELISA: 1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa cells using CDC2/CDC2 (pThr160) antibody.



Immunohistochemical analysis of breast cancer tissue using CDC2/CDC2 (pThr160) antibody.