


Product Datasheet

p62 Dok (Phospho-Tyr398) antibody (orb222983)

Description	Rabbit polyclonal antibody to NUP89
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q99704
Dilution Range	WB: 1:500-1:3000, ELISA: 1:20000
Expiration Date	12 months from date of receipt.



Western blot analysis of K562 cells usin...