



Product Datasheet

eIF4E (Phospho-Ser209) antibody (orb222991)

Description	Rabbit polyclonal antibody to eIF4E
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P06730
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, ELISA: 1:20000
Expiration Date	12 months from date of receipt.



Western blot analysis of 293 cells using...



Immunohistochemical analysis of human br...