

Product Datasheet

Estrogen Receptor-alpha (Phospho-Ser106) antibody (orb222998)

Description	Rabbit polyclonal antibody to Estrogen Receptor-alpha
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P03372
Dilution Range	WB: 1:500-1:3000, ELISA: 1:10000
Expiration Date	12 months from date of receipt.