



Product Datasheet

HER2 antibody (orb223022)

Description	Rabbit polyclonal antibody to HER2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P04626
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, ELISA: 1:10000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of breast c...



Western blot analysis of HeLa cells and ...