



Product Datasheet

c-Jun (Phospho-Ser243) antibody (orb223037)

Description	Rabbit polyclonal antibody to c-Jun
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P05412
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, ELISA: 1:20000
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa cells using c-Jun (Phospho-Ser243) antibody.



Immunohistochemical analysis of human breast tissue using c-Jun (Phospho-Ser243) antibody.