


Product Datasheet

JunB (Phospho-Ser79) antibody (orb223043)

Description	Rabbit polyclonal antibody to JunB
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IP
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P17275
Dilution Range	IHC-P: 1:50-1:100, ELISA: 1:20000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of breast c...