


Product Datasheet

MARCKS (Phospho-Ser163) antibody (orb223049)

Description	Rabbit polyclonal antibody to MARCKS
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P29966
Dilution Range	IHC-P: 1:50-1:100, ELISA: 1:10000
Expiration Date	12 months from date of receipt.



Immunohistochemical
analysis of human
br...