

Product Datasheet

p95 (Phospho-Ser343) antibody (orb223060)

Description	Rabbit polyclonal antibody to p95
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	O60934
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, ELISA: 1:40000
Expiration Date	12 months from date of receipt.