



Product Datasheet

p53 (Phospho-Ser15) antibody (orb223073)

Description	Rabbit polyclonal antibody to p53 (Phospho-Ser15)
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	N-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P04637
Dilution Range	IHC-P: 1:50-1:100, ELISA: 1:10000, WB: 1:500-1:3000
Expiration Date	12 months from date of receipt.



Immunohistochemistry analysis of human b...



Western blot analysis of COS7 cells usin...