

Product Datasheet

p53 (Phospho-Ser33) antibody (orb223074)

Description	Rabbit polyclonal antibody to p53 (Phospho-Ser33)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P04637
Dilution Range	IHC-P: 1:50-1:100, ELISA: 1:10000, WB: 1:500-1:3000
Expiration Date	12 months from date of receipt.