



Product Datasheet

p70 S6 Kinase (Phospho-Ser424) antibody (orb223079)

Description	Rabbit polyclonal antibody to p70 S6 Kinase
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P23443
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:100, ELISA: 1:10000
Expiration Date	12 months from date of receipt.



Western blot analysis of 293 cells using...



Immunohistochemical analysis of human br...