


Product Datasheet

PKR (Phospho-Thr451) antibody (orb223086)

Description	Rabbit polyclonal antibody to PKR
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P19525
Dilution Range	IHC-P: 1:50-1:100, ELISA: 1:20000
Expiration Date	12 months from date of receipt.



Immunohistochemical
analysis of human
co...