

Product Datasheet

S6 Ribosomal Protein (Phospho-Ser235) antibody (orb223097)

Description Rabbit polyclonal antibody to S6 Ribosomal Protein

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications WB

Immunogen C-terminal

Concentration 1mg/ml


Note For research use only

Clonality Polyclonal

Uniprot ID **P62753**

Dilution Range WB: 1:500-1:1000, ELISA: 1:20000

Expiration Date 12 months from date of receipt.



Western
blot
analysis
of 293
cells
using...