

Product Datasheet

Src (Phospho-Tyr418) antibody (orb223102)

Description	Rabbit polyclonal antibody to Src
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P12931
Dilution Range	WB: 1:500-1:1000, ELISA: 1:5000
Expiration Date	12 months from date of receipt.