



Product Datasheet

Stathmin 1 (Phospho-Ser37) antibody (orb223115)

Description	Rabbit polyclonal antibody to Stathmin 1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	N-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P16949
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:100, ELISA: 1:20000
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat cells us...



Immunohistochemical analysis of human br...