



Product Datasheet

VASP (Phospho-Ser157) antibody (orb223130)

Description	Rabbit polyclonal antibody to VASP
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P50552
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, ELISA: 1:4000
Expiration Date	12 months from date of receipt.



Western blot analysis of NIH-3T3 cells u...



Immunohistochemical analysis of human to...