



Product Datasheet

AKT antibody (orb223171)

Description	Rabbit polyclonal antibody to AKT
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9Y243 , P31751 , P31749
Dilution Range	WB: 1:1000-1:3000, IHC-P: 1:50-1:100, ELISA: 1:10000
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa cells using AKT antibody (orb223171).



Immunohistochemistry analysis of human breast tissue using AKT antibody (orb223171).