


Product Datasheet

BRCA1 antibody (orb223173)

Description	Rabbit polyclonal antibody to BRCA1
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P38398
Dilution Range	IHC-P: 1:50-1:100, ELISA: 1:4000
Expiration Date	12 months from date of receipt.



Immunohistochemistry analysis of human b...