



Product Datasheet

EGFR antibody (orb223178)


Description	Rabbit polyclonal antibody to EGFR
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P00533
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, ELISA: 1:20000
Expiration Date	12 months from date of receipt.



Western blot analysis of A549 cells and ...



Immunohistochemistry analysis of human b...



Western blot analysis of HT-29 cells usi...