



Product Datasheet

I kappa B-alpha antibody (orb223182)

Description	Rabbit polyclonal antibody to I kappa B-alpha
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	N-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P25963
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, ELISA: 1:10000
Expiration Date	12 months from date of receipt.



Western blot analysis of MCF7 cells using I kappa B-alpha antibody.



Immunohistochemistry analysis of human breast tissue using I kappa B-alpha antibody.