



Product Datasheet

Tau antibody (orb223201)

Description	Rabbit polyclonal antibody to Tau
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P10636
Dilution Range	WB: 1:500-1:3000, ELISA: 1:20000
Expiration Date	12 months from date of receipt.



Western blot analysis of 293 cells using...



Western blot analysis of HeLa cells usin...