


Product Datasheet

SOX6 antibody (orb223222)

Description	Rabbit polyclonal antibody to SOX6
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P35712
Dilution Range	WB: 1:500-1:3000, ELISA: 1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat cells us...