


## Product Datasheet

### BLCAP antibody (orb223229)

|                            |                                     |
|----------------------------|-------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to BLCAP |
| <b>Species/Host</b>        | Rabbit                              |
| <b>Reactivity</b>          | Human, Mouse, Rat                   |
| <b>Conjugation</b>         | Unconjugated                        |
| <b>Tested Applications</b> | IF                                  |
| <b>Immunogen</b>           | Internal                            |
| <b>Concentration</b>       | 1mg/ml                              |
| <b>Note</b>                | For research use only               |
| <b>Clonality</b>           | Polyclonal                          |
| <b>Uniprot ID</b>          | <b>P62952</b>                       |
| <b>Dilution Range</b>      | IF/ICC: 1:100-1:500, ELISA: 1:40000 |
| <b>Expiration Date</b>     | 12 months from date of receipt.     |



Immunofluorescence  
analysis of 3T3  
cells...