



Product Datasheet

LMTK2 antibody (orb223232)


Description	Rabbit polyclonal antibody to LMTK2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q8IWU2
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, IF/ICC: 1:100-1:500, ELISA: 1:20000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HepG2 cells...



Western blot analysis of LOVO cells using LMTK2 antibody.



Immunohistochemistry analysis of human brain tissue.