



Product Datasheet

ADCY5/6 antibody (orb223243)

Description	Rabbit polyclonal antibody to ADCY5/6
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	O43306, O95622
Dilution Range	IHC-P: 1:50-1:100, IF/ICC: 1:100-1:500, ELISA: 1:20000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...



Immunohistochemistry analysis of human b...